

CATALOG



505-881-1376 WWW.OSSOTANIUM.COM

CONTENTS

Implants

	Bone Level	2
	Tissue Level	3
	One-Piece Abutment	4
	One-Piece O-Ball	5
	OMNI-Tight	6
Restorat	ive Components	
	3.5mm Bone Level	7-8
	4.3 & 5.0mm Bone Level	9
	Bone Level & Tissue Level	10-11
	Tissue Level	12

One-Piece _____13

OMNI-Tight	 14-15

Surgical Placement Tools & Instruments

Surgical Kits	16
Implant & Abutment Drivers	17
Surgical Placement Instruments	18-19
Implant Motors & Handpieces	20
Regenerative Materials	21



QMS CERTIFIED

BIOTANIUM®

Biotanium is a nanostructured variant of pure titanium. Biotanium is a superior material due to undergoing a superplastic deformation process which reduces the grain size 500-1000 times. This nanostructuring process improves the strength of implant and decreases your patients' healing time.

- Pure titanium

- Does not contain aluminum or toxic vanadium
- Surface treatment: acid rinsed; surface textured with passivated MCD apatitic abrasive (300-800 nm) multi-phase calcium phosphate, leaving a residue-free surface after passivation
- Faster osseointegration due to higher surface morphology of grain size
- Stronger than standard alloys
- Decreased bacterial adhesion due to nonporous surface, also decreasing the risk of peri-implantitis

TITANIUM

Pure titanium implants are made of high grade surgical titanium. When you purchase implants and components from Ossotanium, have confidence that you're providing safe, high quality products for your patients.

- Pure titanium
- Does not contain aluminum or toxic vanadium
- Surface treatment: acid-rinsed; surface textured with passivated MCD apatitic abrasive (300-800 nm) multi-phase calcium phosphate, leaving a residue-free surface after passivation

Components marked with an asterisk (*) may require additional shipping charges and time for delivery.

BONE LEVEL IMPLANTS "OSSO-TIGHT"

Compatible restorative components for 3.5mm diameter on pp. 7-8, 4.3mm & 5.0mm diameter on pp. 9-11

Compatible with Screw Retained Abutments

<u>Thread ØLength</u>		<u>Single</u>	<u>Kit</u>	
3.5mm	Х	8.5mm 10mm 11.5mm 13mm 14.5mm 16mm	B3585S B3510S B35115S B3513S B35145S B3516S	B3585K B3510K B35115K B3513K B35145K B3516K
4.3mm	х	8.5mm 10mm 11.5mm 13mm 14.5mm 16mm	B4385S B4310S B43115S B4313S B43145S B4316S	B4385K B4310K B43115K B4313K B43145K B4316K
5.0mm	x	7.5mm 8.5mm 10mm 11.5mm 13mm 14.5mm	B5075S B5085S B5010S B50115S B5013S B50145S	B5075K B5085K B5010K B50115K B5013K B50145K

Drivers (p. 17)	Short	Long
Implant Driver (3.5mm diameter)	PI533	PI534
Implant Driver (4.3mm and wider diameter)		PI536
Healing Screw & Abutment Torx Drivers	PI537*	PI538
Hex O-Ball Driver (for O-ball abutments)	PI541	PI542

* short Torx driver only recommended for shortened abutments in tight spaces

Locator[®] Placement Tool

08393



Implant Single Includes:



TITANIUM \$150 BIOTANIUM \$225

Implant

Implant Kit

Includes: Implant Analog Healing Abutment (2.0mm height) Torx Screw (for abutment) One Abutment of Your Choice (see list below)

TITANIUM \$250



+ One Abutment of Your Choice from the following: Angled Abutment, Straight Abutment, Transfer Abutment, Prepable Abutment, or Ball Abutment (see pp. 7-11)

TISSUE LEVEL IMPLANTS "PERIO-TIGHT"

Compatible restorative components on pp. 10-12

Compatible with Screw-Retained Abutments

<u>Threa</u>	<u>nd</u>	<u>ØLength</u>	<u>Sin</u>	<u>gle</u>	<u>Kit</u>
3.6mm	х	11mm 13mm 15mm	T3611 T3613 T3615	S	T3611K T3613K T3615K
3.7mm	х	11mm 13mm 15mm	T3711 T3713 T3715	3S	T3711K T3713K T3715K
4.1mm	х	11mm 13mm 15mm	T4111 T4113 T4115	S	T4111K T4113K T4115K
4.6mm	х	9mm 11mm 13mm 15mm	T4609 T4611 T4613 T4615	IS BS	T4609K T4611K T4613K T4615K
5.1mm	х	9mm 11mm 13mm 15mm	T5109 T5111 T5113 T5115	IS BS	T5109K T5111K T5113K T5115K
5.8mm	х	8mm 10mm 12mm	T5808 T5810 T5812)S	T5808K T5810K T5812K
Drivers (p.	17)		Short	Long	
Implant D	ivorc	(all diamotors)	DT525	DT526	

Implant Drivers (all diameters)PI535PI536Healing Screw & Abutment Torx DriverPI537*PI538Hex O-Ball Driver (for O-ball abutments) PI541PI542

* short Torx driver only recommended for shortened abutments in tight spaces

Locator[®] Abutment Placement Tool 0

08393



3.7mm



4.1mm and wider



Implant Single

Includes:

Implant Emergence Profiler (2.7mm) Healing Screw (1.2mm height)

TITANIUM \$150

BIOTANIUM \$225

Implant Kit

Includes:

Implant Emergence Profiler (2.7mm) Healing Screw (1.2mm) Analog Torx Screw (for abutment) One Abutment of Your Choice (see list below)

TITANIUM \$250



+ One Abutment of Your Choice from the following: Angled Abutment, Straight Abutment, Transfer Abutment, or Ball Abutment (see pp. 10-12)

ONE-PIECE ABUTMENT IMPLANTS

Compatible restorative components on p. 13

<u>Thread ØLength</u>		<u>Single</u>	<u>Kit</u>	
2.5mm	х	10mm (3.0mm Platform)	A2510S	A2510K
3.0mm	х	11mm 13mm 15mm	A3011S A3013S A3015S	A3011K A3013K A3015K
3.3mm	х	11mm 13mm 15mm	A3311S A3313S A3315S	A3311K A3313K A3315K
3.8mm	х	11mm 13mm 15mm	A3811S A3813S A3815S	A3811K A3813K A3815K
4.3mm	х	11mm 13mm 15mm	A4311S A4313S A4315S	A4311K A4313K A4315K

Drivers (p. 17)	Short	Long
One-Piece Abutment Implant Driver	PI543	PI544



Implant Single



TITANIUM **\$79** BIOTANIUM **\$97**


Implant Kit

Includes: Implant Analog Waxing Sleeve



ONE-PIECE O-BALL IMPLANTS

Compatible restorative components on p. 13



OMNI-TIGHT[®] IMPLANTS

Compatible restorative components on pp. 14-15

Tissue Level

Healing Screw Driver

Compatible with Cementable Abutments

<u>Thread</u>	<u>ØLength</u>	<u>Sir</u>	<u>ngle</u>	<u>Kit</u>
3.5mm x	11mm 13mm 15mm	O35 O35 O35	13S	O3511K O3513K O3515K
4.0mm x	11mm 13mm 15mm	O40 O40 O40	13S	O4011K O4013K O4015K
4.5mm x	9mm 11mm 13mm 15mm	045 045 045 045	11S 13S	O4509K O4511K O4513K O4515K
5.0mm x	9mm 11mm 13mm 15mm	050 050 050 050	11S 13S	O5009K O5011K O5013K O5015K
6.0mm x	8mm 10mm 11mm 12mm 13mm 15mm	060 060 060 060 060 060	10S 11S 12S 13S	O6008K O6010K O6011K O6012K O6013K O6015K
Placement Dr	ivers (p. 17)	Short	Long	
Implant Drive	er (Octagram)	PI545	PI546	
Healing Cap I	nsertion Tool	PI547	PI548	

PI549





Implant Kit Includes:

Implant Emergence Profiler (2.7mm) Healing Cap (1.2mm) Transfer Device

> TITANIUM **\$150** BIOTANIUM **\$225**

Implant Kit with Restorative Kit (see pp. 14 & 15)

Includes:

Implant Emergence Profiler (2.7mm) Healing Cap (1.2mm) Transfer Device 0°-16° Angled Cementable Abutment Analog Waxing Sleeve

> TITANIUM **\$250** BIOTANIUM **\$325**

3.5mm BONE LEVEL RESTORATIVE COMPONENTS

7

Placement of screw retained restorative components on this page require a Torx driver (p. 17) unless otherwise stated. Healing Abutments and Screws

ĊAE req

Components compatible with 3.5mm diameter bone level implants on p. 2

Analog			\$40
H I	<u>Diameter</u>	<u>Color</u>	
	3.5mm	Purple	BL330
Angled Abutmen	t		\$100
screw retained	Angulation	Platform Heig	<u>ght</u>
includes Torx screw	0°-16° 0°-16°	1.5mm 3.5mm	BL333 BL354
Transfer Abutmen	it		\$85
screw retained includes Torx screw	<u>Angulation</u>	<u>Platform Heic</u>	<u>ght</u>
₩.	0°	1.5mm	BL329
Prepable Abutme	nt		\$85
screw retained	Angulation	Platform Heig	<u>ght</u>
includes Torx screw	0°	1.5mm	BL341
Straight Abutmen	it 👔		\$85
screw retained includes Torx screw		<u>Platform Heic</u> 1.5mm 3.5mm	<u>ght</u> BL355 BL356

requires Torx driver (p. 17)	\$45	
	Platform Height	
	0.5mm (screw)	BL332
	2.0mm (abutment)	BL331
	3.5mm (abutment)	BL342
Castable Abutment		\$25
screw retained		
includes Torx screw	Rexolite®	BL348
*Multi-Unit Abutment	ts	\$130

Multi-Unit Abutments

screw retained, compatible drivers found on p. 17

Torx Driver required for Multi-Unit Abutments with angulation 17° and 30° Straight Multi-Unit Driver requried for Multi-Unit Abutments with angulation of 0°

	Angulation	<u>Platform Heig</u>	<u>ht</u>
	0° (straight)	1.5mm 2.5mm 3.5mm	8000875 8000876 8000877
0° angulation	17°	2.5mm 3.5mm	8000880 8000881
17° angulation	30°	3.5mm 4.5mm	8000885 8000886

Multi-Unit Attachments requires Torx Driver (p. 17)		
Rexolite® Multi-Unit Sleeve with Screw	BL349	\$35
Multi-Unit Ball	BL351	\$25

3.5mm BONE LEVEL RESTORATIVE COMPONENTS

Drivers for placement noted in description if needed. Drivers are found on p. 17.

Components compatible with 3.5mm diameter bone level implants on p. 2



* Locator [®] Abu requires Locator [®] Placeme	tments Int Tool (p. 17)		\$119.50	
* Locator [®] Atta requires Locator [®] Placeme	1.0mn 2.0mn 3.0mn 4.0mn 5.0mn 6.0mn	י ז ז ז	02069 02070 02071 02072 02073 02074	
	1 Denture Cap with 1 Blue Replacement 1 Pink Replacement 1 Clear Replacemen 1 Block Out Spacer	Cap (1.5 lbs) Cap (3 lbs)	ing Male	
$\bigcirc \bigcirc \bigcirc \bigcirc$	2 Kit Pack	8519-2	\$63.50	
	10 Kit Pack	8519-10	\$284	
* Engaging UC	LA (requires Torx	Driver, p. 17)		
includes Torx screw	Plastic	8000810	\$50	
	Gold/Plastic	NL-PGA	\$175	
* Non-Engaging UCLA (requires Torx Driver, p. 17)				
includes Torx screw	Plastic	8000815	\$50	
	Gold/Plastic	8000820	\$175	

4.3mm & 5.0mm BONE LEVEL RESTORATIVE COMPONENTS

Placement of screw retained restorative components require a Torx driver (p. 17).

Components compatible with 4.3mm & 5.0mm diameter bone level implants on p. 2



BONE & TISSUE LEVEL RESTORATIVE COMPONENTS

for Bone Level Implants of Diameters 4.3mm and 5.0mm, and 3.7mm and Larger Tissue Level Implants



Bone Level Restorative Kit Includes: Analog. Transfer Abutment, and Torx screw

	,		- ,	Diameter	\$120
Includes:	BL329,	BL330		3.5mm	RC35
	BL319,	BL326		4.3mm	RC43
	BL319,	BL327		5.0mm	RC50
*Engag	ing UC	CLA (r	equires Torx o	driver, p.17)	
Includes Torx	Screw		Plastic	8000910	\$50
Research and the second se		#12	Gold/Plastic	SL-PGA	\$175
*Non-E		ng U(CLA (require Plastic	es Torx Driver, 8000915	, p. 17) \$50
			Gold/Plastic	8000915	\$175
Castab	le O-Ba	all Ex	tender		\$25
	2	<u>Exten</u>	<u>sion Height</u>	Drive	
Acrylic		4.75m	ım	Hex	OP219
Castab	le Abu	tme	nt		\$25
Screw Reta Includes To				Rexolite®	BL339

BONE & TISSUE LEVEL RESTORATIVE COMPONENTS for Bone Level Implants of Diameters 4.3mm and 5.0mm, and 3.7mm and Larger Tissue Level Implants

* Locator [®] Abutments Requires Locator Placement Tool (p. 17)		\$119.50		evel Restorat		\$120
	<u>Platform Height</u>		Includes:	TL428, TL429	3.6mm	RC36
				TL426, TL420	3.7mm	RC37
	1.0mm	02075		TL426, TL421	4.1mm	RC41
				TL426, TL422	4.6mm	RC46
	2.0mm	02076		TL426, TL423	5.1mm	RC51
	3.0mm	02077		TL426, TL424	5.8mm	RC58
	4.0mm	02078	*Multi-I	Jnit Abutmer	nts	\$130
	5.0mm	02079				4130

02080

Screw Retained, drivers found on p. 17

Torx Driver Required for Multi-Unit Abutments with angulation 17° and 30° Straight Multi-Unit Driver required for Multi-Unit Abutments with angulation of 0°

	Angulation	Platform Heig	<u>lht</u>
	0° (straight)	1.5mm 2.5mm 3.5mm	8000975 8000976 8000977
0° angulation		4.5mm	8000978
	17°	2.5mm 3.5mm	8000980 8000981
17° angulation	30°	3.5mm 4.5mm	8000985 8000986

Multi-Unit Attachments

Requires Torx Driver (p. 17)		
Rexolite® Multi-Unit Sleeve with Screw	BL349	\$35
Multi-Unit Ball	BL351	\$25

* Locator Attachment Kit Requires Locator Placement Tool (p. 17)

Kit Includes:	1 Denture Cap with 1 Blue Replacement 1 Pink Replacement 1 Clear Replacemen 1 Block Out Spacer	Cap (1.5 lbs) Cap (3 lbs)	ng Male
	2 Kit Pack	8519-2	\$63.50
	10 Kit Pack	8519-10	\$284
- Ball Abutn	\$80		

6.0mm

Requires Hex O-Ball Driver (p.17) Kit Includes Cap & O-Rings

 \mathbf{O}



Platform Height	<u>Drive</u>	
0.5mm	Hex	BL320
1.5mm	Hex	BL321
2.5mm	Hex	BL322
3.5mm	Hex	BL323
4.5mm	Hex	BL324

TISSUE LEVEL RESTORATIVE COMPONENTS

Components compatible with screw retained 3.5 mm bone are compatible with tissue implants of diameter 3.6mm on pp. 2 & 3 Components compatible with all screw retained tissue level implants on p. 3

Abutments

```
$100
```

Screw Retained, Includes Torx Screw

		0°-16°	1.5mm	(Cementable)	
--	--	--------	-------	--------------	--

Waxing Sleeves

\$20

Compatible with TL400 abutment	Implant Diameter	
	3.7mm 4.1mm 4.6mm 5.1mm	TL405 TL406 TL407 TL408
Height 6 mm	5.8mm	TL400
Transfer Device		\$25



Compatible with the TL400 abutment

TL404


\$45

Healing Screw (requires Torx driver) Platform Height

Ÿ

1.2mm

TL425



Color

Yellow

Purple

Diameter

TL429

TL420

TL421

TL422

TL423

TL424

3.6mm

3.7mm

4.1mm

4.6mm

5.1mm

5.8mm

	Gold
	Green
	Orange
	Red

Emergence Profiler

2.7mm Height Snap onto Healing Screw (TL425)

<u>Diameter</u>	<mark>\$25</mark> Single	\$150 <u>10 Pack</u>
3.7mm	TL410	TL410-10
4.1mm	TL411	TL411-10
4.6mm	TL412	TL412-10
5.1mm	TL413	TL413-10
5.8mm	TL414	TL414-10

ON	E-PIE	ECE RE	STOR	TIVE (COMPONEN	ITS	
Components compatible with One Piece implants O-Ball A on pp. 4 & 5					nalog		\$35
		\$35	12	1.8mm Single Ball	<u>Hex</u> OP223	<u>Square</u> 0P221	
		<u>meter</u>			2.2mm Single Ball	<u>Hex</u> 0P203	
		nm & 3.0mm	OP200		2.2mm Double Ball	Diameter (He	
		nm & 3.8mm	OP201			3.0mm 3.3mm	OP218 OP222
	4.3r	nm	OP202	O-Ball V	Vaxing Sleeve		\$35
Abutment Implan Waxing Sleeve Pro			\$20			<u>Drive</u>	
			32U			Square	0P224
		<u>neter</u>		Square			
		nm & 3.0mm	OP205	O-Ring & Cap Kit			\$35
	3.3n 4.3n	nm & 3.8mm	OP206 OP207	Kit Includes	: 1 Red O-Ring	<u>1.8mm Ball</u>	<u>2.2mm Ball</u>
O-Ball Titanium				0	1 Black O-Ring 2 Red Block Out Spacers	0P214	OP208
Abutment		Drive	\$35	* 0	1 Housing Cap		
		Hex	OP213	0	0		
				O-Ring	Replacements		\$35
O-BallEmergence	Profiler	\$35	\$150	-		<u>1.8mm Ball</u>	2.2mm Ball
	<u>Drive</u>	Single	<u>10 Pack</u>	Pack of 10 O-	Rings Red 🔘	OP215	OP209
	Hex	OP212	OP212-10		Black	OP217	OP211
Hex Square	Square	OP225	OP225-10	Mixed Pack o	of 5 Black and 5 Red O-Rings	OP216	OP210
			1	3			

OMNI-TIGHT® RESTORATIVE COMPONENTS

Components compatible with OMNI-Tight implants on p. 6

Cementable Post & Core Abutments

	0°-16° Angled	<u>Diameter</u>	\$80
B		3.5mm	OT100
T		4.0mm & wider	OT101
	25° Angled		\$100
		3.5mm	OT102
T		4.0mm & wider	OT103
	0°-16° Angled with F	Platform	\$80
		3.5mm	OT104
		4.0mm & wider	OT105
	0°-16° Angled with 2	2.0mm Extension	\$80
		3.5mm	OT108
		4.0mm & wider	OT109
	0°-16° Angled with 2	2.5mm Extension	\$80
		3.5mm	OT106
		4.0mm & wider	OT107

Healing Cap (r	\$40		
-		<u>Diameter</u>	
		3.5mm	0T116
1.2mm Height		4.0mm & wider	OT117
Healing Screw	V (requires heali	ng screw driver, p. 17)	\$40
		<u>Diameter</u>	
Hade Land		3.5mm	OT118
1.2mm Height		4.0mm & wider	OT119
Emergence P	rofiler	\$30	\$150
	<u>Diameter</u>	<u>Single</u>	<u>10 Pack</u>
	3.5mm	OT122	OT122-10
	4.0mm	OT123	OT123-10
2.7mm Height	4.5mm & 5.0m	m 0T124	OT124-10
	6.0mm	OT125	OT125-10
Restorative K	it		\$110
Analog	·		
Waxing	g Sleeve (not com	patible with the 25° abut	ment)
		<u>Diameter</u>	
Ħ		3.5mm	LC35
		4.0mm	LC40
		4.5mm & 5.0mm	LC45-50
		6.0mm	LC60

OMNI-TIGHT® RESTORATIVE COMPONENTS

Components compatible with OMNI-Tight implants on p. 6

Waxing Sleeve		\$20
Waxing sleeve is not comp	patible with the 25° angled a	butment
	<u>Diameter</u>	
	3.5mm	OT130
	4.0mm	0T131
6	4.5mm & 5.0mm	0T132
	6.0mm	OT133
Analog	Diameter	\$40
	Diameter	
Ħ	3.5mm	0T126
	4.0mm	0T127
	4.5mm & 5.0mm	0T128
	6.0mm	OT129
Transfer Device		\$30
507	<u>Diameter</u>	
32.3	3.5mm	OT120
	4.0mm & wider	0T121

O-Ball Ceme	entable Abutm	nent	\$75
H	Dia	ameter	
	3.5r	mm	0T110
2.2mm ball	4.0r	mm & wider	0T111
O-Ball Exter	der		\$25
for cementable O	ball abutment		
X	3.75mm Castable (<i>i</i>	Acrylic)	0T114
	3.75mm Weldable	(Titanium)	OT115
O-Ring & Ca	p Kit		\$35
for cementable ba	all abutment (2.2mm	ball)	
2 Re	ed O-Ring ack O-Ring ed Block Out Spacers ousing Cap		OP208

SURGICAL PLACEMENT KITS

Surgical Placement Kit \$2,500

for Free-Handed Procedures



Includes:

Autoclavable Tray (maximum 125° C) **Tissue Punches:** 3.5mm, 4.0mm, 4.5mm, 5.0mm, 5.5mm, & 6.5mm Set of Drills, each 16mm long: 1.8mm, 2.4mm, 2.9mm, 3.4mm, 3.9mm, 4.4mm, 4.8mm, & 5.7mm Two Torx Drivers (short & long) Two Drivers for 3.5mm Bone Level Implants (short & long) Two Drivers for 3.7mm & Wider Bone and Tissue Level Implants (short & long) Two Drivers for One Piece Abutment Implants (short & long) Two Drivers for Hex O-Ball Implants (short & long) Lindemann Sidecutting Bur Drill Extender Short Motor Latch Adapter Thumb Knob Paralleling Pins (1 pair) **Double-Ended Ratchet Wrench** Torque-Tester Wrench (22/40 N · cm)

Guided Surgical Placement Kit \$3,500

for Free-Handed & Surgically Guided Procedures

PK2

PK3

Includes PK1 Surgical Kit, Plus:

Extra Set of Drills, each 21mm long: 1.8mm, 2.4mm, 2.9mm, 3.4mm, 3.9mm, 4.4mm, 4.8mm, & 5.7mm Incremental Torque Wrench (15-60 N · cm)

Bone Expansion & Sinus Lift Kit \$650

Includes:

Autoclavable Tray 1.8mm Pilot Drill (16mm length) Bone Expander Screws: 3.0mm, 3.5mm, 4.0mm, & 4.5mm Sinus Lift Screws: 3.0mm, 3.5mm, 4.0mm, & 4.5mm Driver Extender Latch Motor Adapter Double-Ended Ratchet Wrench





IMPLANT & ABUTMENT DRIVERS

IMPLANT DRIVERS

Bone & Tissue Level Implant Drivers					
	Implant Diameter	<u>Short</u>	<u>Long</u>		
	3.5 mm (bone level)	PI533	PI534		
	3.7mm & wider	PI535	PI536		

O-Ball Drivers (for O-ball implants & O-ball abutments)

<u>Drive</u>	<u>Short</u>	Long	
Hex	PI541	PI542	\$65
Square	PI539	PI540	\$50

O-Ball Latch Motor Adapter


For O-ball implants & O-ball abutments

	<u>Drive</u> Square Hex	<u>Short</u> PI579 PI577	<u>Long</u> PI580 PI578	\$60 \$60
One-Piece Ab	outmer	nt Impla	nt Drive	r \$70
			<u>Short</u>	Long
			PI543	PI544
OMNI-Tight [®] I	mplan	t Driver	S \$65	5
		<u>Short</u>	Lon	g
	PI54	5 PI5	646	

ABUTMENT DRIVERS


Torx Abutment Driver For healing screws, healing abutments, and	d screw retained abut	\$50 ments
	<u>Short</u> * PI537	<u>Long</u> PI538
* short Torx driver only recommended for	shortened abutments	in tight spaces
*Locator [®] Placement Too	ol	\$82
		08393

OMNI-Tight[®] Healing Cap Insertion Tool






SURGICAL PLACEMENT INSTRUMENTS



Bone Expander & Sinus Lift Screws \$50





*Incremental Torque Wrench				
	PI557			
15-60 N · cm				

SURGICAL PLACEMENT INSTRUMENTS

Pilot Tissue Punch						\$55
			Diame	<u>eter</u>		
			1.8mn	n		PI500
Tissue Punch						\$75
T	Diame	<u>eter</u>	<u>Color</u>			
		3.5mr	n	Blue		PI501
		4.0mr	n	Red		PI502
		4.5mr	n	Black		PI503
		5.0mr	n	Yellov	V	PI504
		5.5mr	n	Green		PI505
U		6.5mr	n	Orang	je	PI506
Drill Gui	des		\$15		\$40	
	<u>Diameter</u>		<u>Single</u>		<u>10 Pac</u>	<u>k</u>
	2.4mm		PI526	5	PI526	5-10
	2.9mm		PI527	7	PI527	7-10
	3.4mm		PI528		PI528-10	
	3.9mm		PI529)	PI529	9-10
	4.4mm		PI530)	PI530	0-10
	4.8mm		PI531	L	PI531	-10
	5.7mm		PI532	2	PI532	2-10

Pilo	t Drill	3		\$65
	1 and 1	<u>Diameter</u> 1.8mm	<u>16mm length</u> PI508	<u>21mm length</u> PI509
Dril	ls			\$115
Ŧ	<u>Diameter</u> 2.4mm 2.9mm	<u>Color</u> Blue Black	<u>16mm length</u> PI510 PI512	<u>21mm length</u> PI511 PI513
	3.4mm 3.9mm	Purple Yellow	PI514 PI516	PI515 PI517
Ţ	4.4mm 4.8mm 5.7mm	Green Orange Red	PI518 PI520 PI522	PI519 PI521 PI523
Lin		Sidecuttir		\$80
]3.	in the second	<u>Diameter</u> 2.1mm	PI507
Drill Extension		11mm Extens	\$45 ion PI525	
		#		

IMPLANT MOTORS & HANDPIECES

To view all available implant motors, hand pieces & warranty information, please visit aseptico.com

\$1,070

*Implant Motors



6000 Series AEU-6000-70V \$3,220 7000 Series AEU-7000E-70V \$4,580

*Powered Trolley Carts

Implant Motors and Handpieces not included



*Impulsion Dental Implant Handpieces

55 N · cm maximum torque compatible with 6000 Series



20:1 Swing Latch Reducer	AHP-85-I	\$990
20:1 Push Button Reducer	AHP-85P-I	\$1,010

*Mont Blanc Dental Implant Handpieces

80 N · cm maximum torque

compatible with 6000 & 7000 Series

or at	20:1 Reducer	AHP-85MB-X	\$1,590
	20:1 Reducer w/ Depth Gauge	AHP-85MB-CX	\$1,650

REGENERATIVE MATERIALS

*NeoMem®

*C-Blast & C-Graft Putty

NeoMem[®] is a resorbable collagen membrane derived from highly purified type I collagen fibers from bovine Achilles tendon.

15mm x 20mm	N1520	\$109
20mm x 30mm	N2030	\$130
30mm x 40mm	N3040	\$191



Demineralized Bone Matrix with Cancellous Bone



<u>Volun</u>	<u>Volume</u>		
0.3cc	CB-P03	\$99	
0.5cc	CB-P05	\$129	
1.0cc	CB-P10	\$180	
3.0cc	CB-P30	\$464	

Demineralized Bone Matrix with CMC Carrier



-	<u>Volume</u>			
	0.3cc	CG-P03	\$99	
	0.5cc	CG-P05	\$129	
	1.0cc	CG-P10	\$180	
_	3.0cc	CG-P30	\$464	





3321 Columbia Dr NE Albuquerque, NM 87107

505-881-1376

www.ossotanium.com REV N Ossotanium products are manufactured by BDI.